
Disposable Electrosurgical Instruments


LS1001, LS1001S/L, LS1002, LS1002S/L, LS1003, LS1003S/L, LS1004, LS1004S/L, LS1011, LS1011S/L, LS1012, LS1012S/L, LS1021, LS1021S/L, LS1022, LS1022S/L, LS1023, LS1023S/L, LS1024, LS1024S/L, LS1025, LS1025S/L, LS1026, LS1026S/L, LS1027, LS1027S/L, LS1028, LS1028S/L, LS1029, LS1029S/L, LS1030, LS1030S/L, LS1031, LS1031S/L, LS1032, LS1032S/L, LS1033, LS1033S/L, LS1034, LS1034S/L, LS1035, LS1035S/L, LS1036, LS1036S/L, LS1038, LS1038S/L, LS1041, LS1041S/L, LS1042, LS1042S/L, LS1043, LS1043S/L, LS1044, LS1044S/L, LS1045, LS1045S/L, LS6001, LS6001S/L, LS6011, LS6011S/L, LS6021, LS6021S/L, LS6022, LS6022S/L, LS6023, LS6023S/L, LS6024, LS6024S/L, LS6025, LS6025S/L, LS6026, LS6026S/L, LS6027, LS6027S/L, LS6028, LS6028S/L, LS6029, LS6029S/L, LS6030, LS6030S/L, LS6031, LS6031S/L, LS6032, LS6032S/L, LS6033, LS6033S/L, LS6034, LS6034S/L, LS6035, LS6035S/L, LS6036, LS6036S/L, LS2002, LS2002S/L, LS2003, LS2003S/L, LS2004, LS2004S/L, LS2202, LS2202S/L, LS2203, LS2203S/L, LS2204, LS2204S/L

Instructions for use


Monopolar Instrument


Bipolar Instrument


Articulating Instrument


5

6

7

8

EN - ENGLISH

BEFORE USING THIS PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

IMPORTANT

- This booklet is designed to assist using this product. It is not a reference for surgical techniques.
- This device was designed, tested and manufactured for single patient use.
- Reuse of this device may lead to its failure and subsequent patient injury.
- Reprocessing and/or re-sterilization of this device may create risk of contamination and patient infection.
- Do not reuse, reprocess or re-sterilize this device.
- Rx only.

INTENDED USE

The device is intended for the coagulation, cutting, dissecting and grasping tissues during a laparoscopic procedure.

CONTRAINdications

- The electrosurgical instruments are not intended for contraceptive coagulation of fallopian tissue, but may be used to achieve hemostasis following transaction of fallopian tube.
- Patient with severe abdominal inflammation, adhesions and other clinically not justified for laparoscopic procedures are contraindicated for the use of this instrument.
- The use of VHMED laparoscopic products is contraindicated whenever laparoscopic surgical techniques are contraindicated for any reason.

DESCRIPTION

- Rated voltage of monopolar surgical instrument: 2000V. Rated voltage of bipolar surgical instrument: 800V;
- Not made with natural rubber latex.

- Shelf life: 3 years.

WARNINGS AND PRECAUTIONS

- Verify the compatibility of parts from different manufacturers prior to conducting the procedure. And make sure the repeat peak voltage not exceed the rated voltage of the instrument.
- Endoscopic surgery should be performed only by physicians who are thoroughly trained in endoscopic techniques and failure modes, precautions and corrective actions in the event of failure.
- A thorough understanding of techniques and principles involved in electrosurgical procedures are necessary to avoid burn and shock hazards to patient and operator.
- The medical literature should be consulted regarding the techniques, contraindications, hazards and complications prior to the performance of these procedures.
- During high frequency operation if operated improperly may cause unintended damage of tissue, gas embolism, low frequency leakage, influence other electronic equipment, etc. Follow all instructions for use required by the generator manufacturer. Consult medical literature or country specific regulations for specific techniques, complications and hazards prior to procedure.
- Do not use a damaged or faulty instrument. Check the instrument for bent, broken, cracked, worn, or separated parts prior to use.
- Ensure device grasping or cutting surfaces are fully visible prior to engaging the electrical current to avoid unwanted damage. Keep the working-end under full and unobstructed visualization during use.
- The instrument should be used cooperated with endoscope and the accessories in the laparoscopic operation.
- This instrument is not water-proof. Do not use this system in any place where water leakage may occur. If any water is sprayed on or into the instrument handle, electric shock may result. If water is accidentally sprayed on or into the instrument, contact VHMED customer service representative.
- Do not use excessive force or in a manner not consistent with normal instrumentation use.
- This instrument can only clamp or cut part of soft tissues. Check integrity of the instruments after each use to avoid left damaged parts in the patient.
- The instrument is not intended for use when endoscopic techniques are contraindicated.
- Disposable devices are intended for single patient use only-DO NOT RESTERILIZE. Dispose of all used or damaged products using standard hospital practices for biohazard control. Instructions for sterilization of re-usable devices are noted in the appropriate sections.
- There is no contraindication directly linked with the instrument till now. Doctor should have the idea whether the instrument can be used on the patient according to the patients' particular situation. For further information about techniques, complications and hazards, please consult the latest medical literature.

POTENTIAL COMPLICATIONS FOR ALL VHMED PRODUCTS

Potential complications associated with the use of laparoscopic devices include but are not limited to organ or vessel damage / perforation, hemorrhage, hematoma, trocar site hernia, and septicemia / infection.

GENERAL INSTRUCTIONS FOR USE (Monopolar Instrument and Bipolar Instrument)

1. Choose the instrument with appropriate tip characteristics according to the requirements of the surgery. If the connection to highfrequency electrosurgical generator is needed, ensure that the product high frequency function is intact. When attaching a high frequency electrosurgical generator to the electrosurgical post located on the top of the instrument, follow the generator manufacturer's instructions for use and setup. **Refer to picture 1**
2. Close the tip and insert the instrument into the patient cavity through trocar. **Refer to picture 2**
3. Rotate the knob to turn the instrument tip to the desired location. **Refer to picture 3**
4. Position tissue inside the blades / jaw tips and squeeze handle to cut it, or position tissue inside the jaws and

squeeze handle to desired clamping tightness. Different closing angle have different locking forces. Choose the correct closing angle. **Refer to picture 4.0**

For instruments with ratchet, choose the right closing angle and press the ratchet in front of the handle to activate or deactivate the ratchet. Please follow the instruction to operate. **Refer to picture 4.1**

GENERAL INSTRUCTIONS FOR USE (Articulating Instrument)

1. Choose the instrument with appropriate tip characteristics according to the requirements of the surgery. **Refer to picture 1**
2. Rotate the lipstick style articulation wheel to retract the instrument tip straight inside the instrument shaft. Select whether connect the instrument to the high frequency electrosurgical generator. If high frequency is required, ensure that the product high frequency function is intact. When connecting a high frequency electrosurgical generator to the electrosurgical post located on the top of the instrument, follow the generator manufacturer's instructions for use and setup. **Refer to picture 2**
3. Insert the instrument into an appropriate sized trocar port. **Refer to picture 3**
4. To rotate the lipstick style articulation wheel, articulate the tip to the desired location. **Refer to picture 4**
5. Adjust the color rotation knob, turn the tip to desired location, lock the rotation knob. **Refer to picture 5**
6. Position tissue inside the jaws and squeeze handle to desired clamping tightness. Different closing angle has different locking forces. Choose the correct angle. Please follow the instruction to operate. For instruments with ratchet, choose the correct angle and press the ratchet in front of the handle to close. Please follow the instruction to operate. **Refer to picture 6**
7. Rotate the lipstick style articulation wheel to retract the tip to the straight position. **Refer to picture 7**
8. Remove the instrument from the trocar port. **Refer to picture 8**

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VH MED (Nantong) Co. Ltd.